

tobiidynavox

Sono Key

Intuitive Communication and Access
to all functions of your AAC device



Power to be You

The key to all functions of your communication device

Intuitive communication and access to all functions of your communication device. Ideal for all communication needs of literate people.



Easy access to your device

Sono Key is a user interface for Tobii Dynavox Communicator that provides streamlined, comprehensive and easy access to your communication device.

It is designed to get literate users up and running quickly and provides complete access to all of the functionality inside your computer and communication device. The intuitive, yet powerful design allows for a wide range of input methods, giving anyone the chance to use modern information technology to enrich their lives.

Control speech and other functions

With Sono Key you can create your own speech-synthesized messages as well as store and print documents. You can also control other key computer functions like calendars, e-mail, chat, SMS, media players, photo libraries and Windows controls.

Customisable in many ways

Personalising input methods is accomplished in a few seconds and you can even choose from several different keyboard layouts including, abc, QWERTY and keyboards optimised for eye control.

Sono Key offers a variety of selection methods including, direct selection by click or touch, selection with adjusted dwell time, scanning and even eye control.

All settings can easily be adjusted with no programming required.

Sono Key features

- Text-to-speech – type and let your computer speak for you.
- Manage text – open and save documents or send them to a printer.
- Quick phrases – store as many premade sentences as you want, sort them into groups, and communicate faster.
- E-mail control – send and receive e-mail with the support of synthetic speech.
- Instant messaging – exchange messages to keep contact with family and friends
- Calendar functionality – browse your calendar, add appointments, and create alerts for important events.
- Windows on-screen keyboard – gives extended functionality and easy access for web browsing.
- Quick Talk – speak without having to change applications.
- Music Player – listen to your favourite music whenever you want.
- Photo Browser – view photos and create your own slideshows.
- Calculator – math made easy.
- Mobile phone functionality - send and

receive SMS messages and make regular phone calls.

- Volume control – quickly increase or decrease volume.
- Environmental control – program environmental control functions for TV, lights, the household and more.
- Alerts – program attention grabbing alert sounds on your computer.
- Screen functions – easily open, close, enlarge or reduce the size of applications.
- Oversized buttons – large and easy to access buttons.

System requirements

Sono Key requires Tobii Dynavox Communicator Version 4.2 Premium or higher.

Some of the functions require a mobile phone subscription and an internet connection.

More information

To find out more about Sono Key, please visit www.tobiidynavox.com or any of our other websites.

© TobiiDynavox™. Illustrations and specifications do not necessarily apply to products and services offered in each local market. Technical specifications are subject to change without prior notice. All other trademarks are the property of their respective owners.

We believe that everyone has the right to live their life to the fullest. Wanting to help people do just that, we provide them with the tools to communicate, become independent and do what they once were able to do or never thought possible.