

Feature group	Feature	Microsoft Eye Control	Windows Control	Gaze Point	Classic Tobii Gaze Interaction Software	Gaze Interaction + Communication software ¹
Keyboard	Gaze controlled typing	✓	✓	✗	✓	✓
	Word prediction	✓	✓	✗	✓	✓
	Text-to-Speech	✓	✗	✗	✗	✓
Mouse control	Moving cursor ²	✗	✓	✓	✓	✓
	Placing cursor	✓	✓	✗	✓	✓
	Zoom-clicking	✗	✓	✗	✓	✓
	Left click	✓	✓	✓	✓	✓
	Gaze Right click, Double Clicktyping	✓	✓	✗	✓	✓
	Reposition the taskbar	✓	✓	N/A	✓	✓
	Scrolling	✗	✓	✗	✓	✓
	Drag & Drop	✗	✓	✗	✓	✓
	Shortcuts ³	✗	✓	✗	✓	✓
	Interaction	Dwell-Free Typing	✗	✓	✓	✓
	Shape Writing	✓	✗	✗	✗	✗
	Switch	✗	✓	✓	✓	✓
	Blink	✗	✗	✗	✓	✓
Calibration	Custom number of calibration points	✗ ⁴	✓	✓	✓	✓
	Image as calibration stimuli	✗	✓	✓	✓	✓
	Video as calibration stimuli	✗	✗ ⁵	✗ ⁵	✓	✓
	Step by step calibration	✗	✓	✓	✓	✓
Other	Region interaction	✗	✓	✓	✓	✓
	Save Calibration Profiles	✗	✓	✓	✓	✓ ⁶
	Gaze enabled settings	✗	✓	✓	✗	✓ ⁶
	Tracking one eye only	✗	✓	✓	✓	✓
	Tutorial / Configuration Guide	✗	✓	✓	✗	✓
	Customizable Taskbar	✗	✓	✗	✗	✓

1. A combination of Tobii Dynavox Eye Gaze Software and Communication Software (Communicator 4, Communicator 5, Compass, Snap + Core First)

2. To play games like Look2Learn, EyeFx, etc. that are controlled by mouse movements

3. Direct and quick access to Windows most used functions

4. 7 pts. through TET Core

5. Coming soon

6. Communicator 5