







# Recommended Tobii Eye Trackers for Microsoft Eye Control

<h3>Tobii Eye Tracker 4C</h3>  <p>PC gaming device that provides simultaneous head &amp; eye tracking functionality</p>	<h3>PCEye Mini</h3>  <p>Navigate &amp; control your computer using only your eyes</p>	<h3>PCEye Plus</h3>  <p>Eye tracking, dictation, switch, &amp; IR in one product</p>	<h3>EyeMobile Plus</h3>  <p>Purpose-built bracket solution for mobile communication &amp; complete hands-free access</p>	<h3>I-Series+</h3>  <p>Speech generating device for effective communication through eye tracking</p>
--	--	---	---	---

## Designed for >

				
---	---	---	---	---

## Eye Tracking Performance >


**97%**  
 Eye trackability of population (world wide)

## Included Gaze Software >

<ul style="list-style-type: none"> <li>none</li> </ul>	<ul style="list-style-type: none"> <li>Windows Control</li> </ul>	<ul style="list-style-type: none"> <li>Windows Control &amp; Virtual Remote for the EyeR Dongle</li> </ul>	<ul style="list-style-type: none"> <li>Windows Control &amp; Virtual Remote</li> </ul>	<ul style="list-style-type: none"> <li>TGIS-Classic Windows Control</li> </ul>
--	---	--	--	--

## Other Compatible Gaze Software >

<ul style="list-style-type: none"> <li>Microsoft Eye Control</li> <li>Gaze Point (free)</li> <li>Windows Control (not included)</li> </ul>	<ul style="list-style-type: none"> <li>Microsoft Eye Control</li> <li>Gaze Point (free)</li> <li>Gaze Viewer 1.1.4 (not included)</li> </ul>	<ul style="list-style-type: none"> <li>Microsoft Eye Control</li> <li>Gaze Point (free)</li> <li>Gaze Viewer 1.1.4 (not included)</li> </ul>	<ul style="list-style-type: none"> <li>Microsoft Eye Control</li> <li>Gaze Point (free)</li> <li>Gaze Viewer 1.1.4 (not included)</li> </ul>	<ul style="list-style-type: none"> <li>Microsoft Eye Control (coming soon)</li> <li>Gaze Viewer 1.1.4 (not included)</li> </ul>
--	--	--	--	---

## Compatible AAC Software >

<ul style="list-style-type: none"> <li>Communicator 5.4 (when used with Windows Control)</li> </ul>	<ul style="list-style-type: none"> <li>Communicator 5.4</li> </ul>	<ul style="list-style-type: none"> <li>Communicator 5.4</li> </ul>	<ul style="list-style-type: none"> <li>Communicator 5.4 (included in AAC bundle version)</li> </ul>	<ul style="list-style-type: none"> <li>Communicator</li> <li>Snap + Core First</li> <li>Compass</li> <li>ALL</li> </ul>
---	--	--	---	---

## Environmental Control (IR), Noise Cancelling Microphone (Dictation), Switch Port build in >

IR	Dictation	Switch	IR	IR	Dictation	Switch
X	X	X	✓	✓	X	✓

## Designed for Assistive Technology >

X	✓	✓	✓	✓
---	---	---	---	---

## Support >

<ul style="list-style-type: none"> <li>Gaming and Developer support via FAQ, email forms, social media and chat</li> </ul>	<ul style="list-style-type: none"> <li>Tobii Dynavox support via our support web page, social media, email, chat (in US), and phone (if purchased or via reseller)</li> </ul>	<ul style="list-style-type: none"> <li>Tobii Dynavox support via our support web page, social media, email, chat (in US), and phone (if purchased or via reseller)</li> </ul>	<ul style="list-style-type: none"> <li>Tobii Dynavox support via our support web page, social media, email, chat (in US), and phone (if purchased or via reseller)</li> </ul>	<ul style="list-style-type: none"> <li>Premium support with Support360<sup>SM</sup> or TD Care (Repair Warranty, Damage Protection and much more)</li> </ul>
--	---	---	---	--